

### Recommendations for Good Lab Results

- Follow directions exactly.
- Know why directions are being followed -that is, the reason for each step of the laboratory procedure.
- Keep a laboratory notebook in which all the directions, results, and calculations are carefully recorded.
- Be well organized in your thinking, in your notebook, and in the laboratory operations.
- Know how to use each item of lab equipment and treat it with great care.
- Know and observe the proper steps for cleaning lab glassware.
- Follow the safety rules.
- Know which laboratory chemicals are hazardous; use them in accordance with safety regulations, and discard them in accordance with hazardous waste regulations.
- Be calm and deliberate in your work.
- Be considerate of others in the laboratory in all that you do (i.e. clean up all messes immediately, leave all equipment the way you would expect to find it, etc ...).

### Estimated Precision of Common Laboratory Apparatus--

Buret: ±0.02 mL for each reading

Analytical Balance: ±0.0001 g

Volumetric Pipets (Class A)	
Volume	Tolerance
5 ml	± 0.1 mL
10 mL	±0.02 mL
15 mL	±0.03 mL
20 mL	±0.03 mL
25 mL	±0.03 mL
50 mL	±0.05 mL
100 mL	±0.08 mL

Volumetric Flasks (TC, Class A)	
Volume	Tolerance
25 mL	±0.03 mL
50 mL	±0.05 mL
100 mL	±0.08 mL
250 mL	±0.12 mL
500 mL	±0.15 mL
1000 mL	±0.30 mL
2000 mL	±0.50 mL

### Grades of Chemicals

Technical or Commercial: indeterminate quality of purity; may be used in preparation of cleaning solution only.

C.P. (Chemically Pure): more refined than technical, but still unknown quality.

USP: Meets minimum purity standards of tolerance set by the United States Pharmacopoeia for contaminants dangerous to human health.

A.C.S. Reagent: High purity; conforms to minimum specifications set by the Reagent Chemicals Committee of the American Chemical Society.

Primary Standard: Highest Purity; required for accurate volumetric analysis (for standard solutions).

Reference Standards: Materials of known composition used most frequently for validation of methods of analysis.

### Concentration of Commercial Reagent-Grade Acids and Bases

Reagent	F.W.	Molarity	% (w/w)	Density (g/mL) @20°C
H <sub>2</sub> SO <sub>4</sub>	98.08	17.6	94.0	1.831
HClO <sub>4</sub>	100.5	11.6	70.0	1.668
HCl	36.46	12.4	38.0	1.188
HNO <sub>3</sub>	63.01	15.4	69.0	1.409
H <sub>3</sub> PO <sub>4</sub>	98.00	14.7	85.0	1.689
HC <sub>2</sub> H <sub>3</sub> O <sub>2</sub>	60.05	17.4	99.5	1.051
NH <sub>3</sub>	17.03	14.8	28.0	0.898

### Rules for Handling Reagents-

- Use the best available grade for analytical work.
- Replace the top of the container immediately after removal or reagent
- Hold stoppers for reagent bottles between your fingers and avoid setting them on the benchtop.
- Never, ever return excess reagent to the bottle.
- Never insert spatulas into a bottle (pour the reagent with a rolling motion from the bottle). There will be occasional exceptions to this rule when reagents "clump".
- Keep your working area neat and dean. Spills are to be wiped up immediately!
- Label EVERYTHING!

### Dixon's Q-Test for Outliers

$$Q_{\text{calculated}} = \frac{|\text{suspect value} - \text{closest value}|}{|\text{highest value} - \text{lowest value}|} = \frac{\text{gap}}{\text{range}}$$

If  $Q_{\text{calculated}} > Q_{\text{critical}}$ , the suspect value should be discarded.

Number of Observations	Q <sub>critical</sub> At 90% confidence
3	0.94
4	0.76
5	0.64
6	0.56
7	0.51
8	0.47
9	0.44
10	0.41

### Grubbs Test for Outliers

$$G_{\text{calculated}} = \frac{|\text{suspect value} - \bar{x}|}{s}$$

$\bar{x}$  = mean

s = sample standard deviation

If  $G_{\text{calculated}} > G_{\text{critical}}$ , the suspect value should be discarded.

Number of Observations	G <sub>critical</sub> At 95% confidence
4	1.463
5	1.672
6	1.822
7	1.938
8	2.032
9	2.110
10	2.176

### Confidence Intervals

$$\mu = \bar{x} \pm \frac{ts}{\sqrt{n}}$$

$\bar{x}$  = mean

s = sample standard deviation

n = number of data points

t = "Student's t"

Degrees of Freedom	t at 95% confidence level	Degrees of Freedom	t at 95% confidence level
1	12.7	7	2.36
2	4.30	8	2.31
3	3.18	9	2.26
4	2.78	10	2.23
5	2.57	15	2.13
6	2.45	∞	1.96